


RBC Kontroll-DNA Panel



- ⇒ Zu Validierungszwecken
- ⇒ Gebrauchsfertig
- ⇒ Vortypisiert für Blutgruppen und HPA
- ⇒ 6 DNAs à 50 µl
- ⇒ Konzentration 40 ng/µl
- ⇒ Artikelnummer 001 K00 006

	F44	F54	F34	12458	12456	12457
AB0	O1O1	O1O1	BB	O1A2	O1O1	O1B
RHD	D	D	D	D	dd	D
RHCE	Cce	Cwce	Cce	Cwe	ce	Cce
Zygotie	Dd	Dd	Dd	DD	dd	Dd
weak D	-	-	type 3	-	-	-
KEL	KEL1	KEL2	KEL2	KEL2	KEL2	KEL2
JK	JK*A	JK*B	JK*A,B	JK*A, B	JK*A, B	JK*A, B
FY	FY*A, B	FY*A, B	FY*A,B	FY*B	FY*A,B	FY*A
MNS	MSs	MNSs	Ns	MSs	MS	MSs
Kp	bb	bb	bb	bb	bb	bb
Lutheran	bb	bb	bb	bb	ab	bb
Diego	bb	bb	bb	bb	bb	bb
Wright	bb	bb	bb	bb	bb	bb
Cartwright	aa	aa	aa	ab	aa	aa
Colton	aa	aa	aa	ab	aa	aa
Knops	aa	aa	aa	aa	aa	ab
Dombrock	bb	bb	bb	ab	bb	bb
HPA 1	aa	aa	aa	aa	aa	aa
HPA 2	aa	aa	aa	aa	aa	aa
HPA 3	ab	ab	aa	ab	ab	ab
HPA 4	aa	aa	aa	aa	aa	aa
HPA 5	aa	aa	aa	aa	aa	aa
HPA 6	aa	aa	aa	aa	aa	aa
HPA 9	aa	aa	aa	aa	aa	aa
HPA 15	aa	aa	bb	ab	ab	ab