



HLA-FluoGene^{NX} ABCDRDQ

Combining 9 HLA loci in one 96-well assay

- Fastest automated platform with the lowest manual hands-on-time
- HLA-A, -B, -C, -DRB1, 3, 4 & 5, DQA1 and DQB1 typing in 1 hour
- Less DNA consumption
- Only 1 manual pipetting step
- Only 59 minutes PCR-Protocol
- Validated for FluoQube real time instrument:
 - Fluorescence signals are detected after each cycle
 - Chemistry and instruments are offered from one hand



ARTICLE NO.	PRODUCT	PLATES/TEST	TESTS/KIT
002 082 010	HLA-FluoGene ^{NX} ABCDRDQ TaqMan typing of HLA-A, -B, -C, -DRB1, 3, 4 & 5, -DQA1 and -DQB1 alleles	1	10

